

# qPCR<sup>BIO</sup> cDNA Synthesis Kit

- Unbiased cDNA synthesis for real-time PCR
- Thermostable reverse transcriptase
- 30 minute protocol

## Features

- Unbiased representation of 5' and 3' mRNA transcript ends
- Sensitive detection of low copy number transcripts
- High cDNA yields from as little as 4pg total RNA
- Simple 2 tube system
- 5x buffer contains anchored oligo(dT), random hexamers, enhancers, dNTPs and MgCl<sub>2</sub>
- 20x thermostable reverse transcriptase blended with RNase inhibitor

## Applications

- cDNA synthesis for real-time PCR analysis
- Low copy number transcripts
- Viral RNA targets
- miRNA targets
- Efficient synthesis from total RNA or poly(A)+ RNA




Figure 1. Broad reverse transcription dynamic range

qPCR<sup>BIO</sup> cDNA Synthesis Kit was used for cDNA synthesis using a 10 fold serial dilution of mouse total RNA from 40pg to 400ng. qPCR was performed using qPCR<sup>BIO</sup> SyGreen Mix amplifying a 122bp fragment of the mouse ACTG gene. Efficiency was measured at 96% across the range tested.

The results demonstrate that qPCR<sup>BIO</sup> cDNA Synthesis Kit efficiently reverse transcribes RNA across a broad dynamic range of substrate.

The qPCR BIO cDNA Synthesis Kit uses the latest developments in reverse transcriptase technology and buffer chemistry to enhance cDNA synthesis speed, yield and representation. The reverse transcriptase, buffer system and combination of random hexamers with anchored oligo(dT) allow for unbiased, efficient and sensitive cDNA synthesis.

High quality cDNA synthesis for downstream qPCR analysis is essential for successful expression studies. Many factors affect cDNA synthesis including the reverse transcriptase, buffer systems, enhancers and priming strategy. The PCR BIO cDNA synthesis mix removes the need for user optimisation of these critical factors.

The modified MMLV reverse transcriptase (RTase) is both thermostable and extremely active. The enzyme is blended with RNase inhibitor preventing degradation of RNA by contaminating RNase. The RTase is not inhibited by ribosomal and transfer RNAs, total RNA is an ideal substrate. The 5x cDNA synthesis mix can be used with up to 0.4µg total RNA.

The relative concentrations of random hexamers and anchored oligo(dT) have been optimised for the generation of cDNA for use in real-time PCR experiments.




Figure 2. Unbiased representation of mRNA ends

a) qPCR BIO cDNA Synthesis Kit was used to synthesise cDNA from mouse liver total RNA. 2 primer pairs were designed against the 5' (red traces) and the 3' (black traces) ends of the 4.2kb mouse CANX transcript. qPCR BIO SyGreen Mix was used for analysis. The primer pairs were 4kb apart and did not show any reverse transcription bias, hence the amplification traces overlap.

b) 2 primer pairs against the 5' (red) and 3' (black) traces of RNS18 gene (1.8kb). Again, no reverse transcription bias was evident.




Figure 3. Thermostable enzyme for high GC%

a) qPCR BIO cDNA synthesis kit was used to synthesise cDNA from mouse liver total RNA at 42°C (red) and also at 55°C (black). A primer pair was designed against GJB2, generating an 84% GC amplicon. qPCR BIO SyGreen Mix was used for analysis. The higher temperature incubation generated more GC rich cDNA than the low temperature incubation.

b) A control amplicon of 55% GC from GAPDH was amplified from the 2 cDNAs described above. For this GC% no advantage of higher temperature incubation was achieved. The yield was slightly lower with the higher temperature.

Catalogue Number	Product Name	Pack size	Presentation
17-700	qPCR BIO cDNA Synthesis Kit	25 Reactions	[1 x 0.1ml] & [1 x 0.025ml]
17-700B		100 Reactions	[4 x 0.1ml] & [1 x 0.1ml]

Distributed by:



Genesee Scientific | 900 Vernon Way | El Cajon | CA 92020  
T: 800-789-5550 | E: support@geneseesci.com  
www.geneseesci.com