

INVICTUS

DROSOPHILA INCUBATOR



INVICTUS Drosophila Incubator Specifications		Cat # 59-400	Cat # 59-400L	Cat # 59-400LH
	Temperature Range	5 - 60°C	5 - 60°C	12 - 27°C
	Temperature Uniformity	± 1°C	± 1°C	± 1°C
	Mite Cycle <i>(60°C decontamination cycle)</i>	Yes	Yes	No
	High/Low Temperature Failsafe	Yes		
	Programmable Lighting Control	No	Yes	Yes
	Additive Humidity Control	No	No	Yes
	Shelves Included <i>(epoxy coated)</i>	6 Shelves		
	Shelf Dimensions	24 x 24in. Holds 4 trays of narrow vials and 3 trays of wide vials/bottles per shelf.		
	E-Coat Corrosion Resistant Coils	Yes		
	Smart Compressor <i>(cat# 59-400LH requires continuously running compressor for optimal humidity control)</i>	Yes	Yes	No
	Door <i>(swinging door, right side hinge)</i>	Glass	Solid	Solid
	Interior Temperature & Humidity Display	Yes		
	Floor Drain Required	No		
	Exterior Color	INVICTUS Royal Blue		
	Exterior Dimensions <i>(W x D x H)</i>	27.9 x 33 x 79in.		
Interior Dimensions <i>(W x D x H)</i>	24 x 25 x 55in.			
Packaging Dimensions <i>(W x D x H)</i>	36 x 43 x 84in.			
L:5:1 Warranty <i>(lifetime on coils, 5 year compressor motor, 1 year parts and labor)</i>	Yes. BEST WARRANTY ON THE MARKET.			

Optional Accessories for INVICTUS	Cat #	Description
	59-401	Extra Shelf, 24 x 24in.
	59-402	Wall Port for Sensor
	59-403	Wall Port for Electrical Plugs
	59-404	Alarm with Relay
	59-405	International Voltage Conversion
	59-406	Continuously Running Refrigeration System - Reduces Ambient RH Levels

White Glove Delivery Service


For an additional charge, Genesee Scientific can arrange for your INVICTUS to be delivered to your lab, installed and for all packing materials to be disposed off. Customer must specify if 'White Glove' delivery service is required at time of order. Genesee Scientific will not be responsible for any redelivery charges.



Cat # 59-400L
& 59-400LH




Cat # 59-400



Interior temperature and relative humidity display. Standard on all INVICTUS models.



INVICTUS Controller



High-Performance Coated Coils



Toll Free: 800.789.5550 Fax: 888.789.0444 Web: www.flystuff.com