

Power Your PCR. Prove Your Results.

Buy a qualifying QTOWER Real Time qPCR Thermal Cycler and receive **2 WEEKS OF PREMIUM CONSUMABLES FREE.**

Get everything you need to thoroughly test your new system with premium reagents and plates — optimized to perform together.



Trusted Performance

QTOWER Thermal Cyclers deliver unmatched precision, reliability, and real-time data you can trust.



Test with Confidence

Receive two weeks of high-quality consumables to validate your workflow from day one.



Quality that Pays Off

Premium consumables. Better data. Stronger results. Fewer repeats.



EXPLORE QTOWER
SYSTEMS

STEP 1: Buy Any QTOWER

QTOWER Real-Time qPCR Thermal Cyclers
High-performance real-time qPCR systems engineered for precision, flexibility, and reproducible results. Optimized for validated workflows and dependable performance across demanding applications.

- ✓ 6-channel multiplexing capability
- ✓ Real-time amplification analysis
- ✓ Compatible with 96- and 384-well workflows

STEP 2: Choose your free bundle

Apply the matching promo code at checkout to receive your preferred qPCR workflow bundle free with your qualifying QTOWER purchase.

Apply Promo Code THERMAL1

SYBR + 96-Well Workflow



- 4x Apex GREEN qPCR Master Mix, 2X, High ROX (42-116PG)
- 2x Olympus 96-Well PCR Plates, Full-Skirted (24-302W)

Apply Promo Code THERMAL2

SYBR + 384-Well Workflow



- 4x Apex GREEN qPCR Master Mix, 2X, High ROX (42-116PG)
- 2x Olympus 384-Well PCR Plates, One Notch (24-305W)

Apply Promo Code THERMAL3

Low-ROX / No-ROX + 96-Well Workflow



- 4x Apex GREEN qPCR Master Mix, 2X, High ROX (42-120PG)
- 2x Olympus 96-Well PCR Plates, Full-Skirted (24-302W)

Apply Promo Code THERMAL4

Low-ROX / No-ROX + 384-Well Workflow



- 4x Apex GREEN qPCR Master Mix, 2X, High ROX (42-120PG)
- 2x Olympus 384-Well PCR Plates, One Notch (24-305W)