

Free Mini Centrifuge with Select Benchmark Homogenizers

Purchase a qualifying
**Benchmark
Homogenizer**
get a Poseidon® Mini
Centrifuge
FREE.

USE PROMO CODE: **GENS26**



HOW IT WORKS:



Add any qualifying
Benchmark Homogenizer
to your cart



Apply promo code
GENS26 at checkout.



Your Poseidon® Mini
Centrifuge will be added
automatically for **FREE!**

YOUR FREE GIFT:

Poseidon® Mini Centrifuge
Compact. Reliable. Ready for Every Lab.



Qualifying Homogenizers include:
31-212, 31-216, 31-214 and 31-214R

Offer Valid Through
September 30, 20206.