

Supplier of the Month Deal

March Special Promo

Promo: **GENMAR**

Valid Through 04.30.2025

Pricing, availability, and specifications subject to change without notice.

MC-24™ Series Microcentrifuges

- Includes COMBI-Rotor™ for tubes and strips
- High speed, up to 16,800xg
- Available with refrigeration
-20°C to + 40°C temp. range



**Touch
Control**



**Quiet
Operation**



**High
Speed**



31-511R
(refrigerated)



31-511
(ambient)

Special Offer:

FREE myFuge Mini Centrifuge
& Mini Vortex Mixer

Receive a 31-500G & 31-503G
with purchase of 31-511 & 31-511R



Item No.	Description	Price
31-511	MC-24 Microcentrifuge w/ COMBI-rotor	\$2,693
31-511R	MC-24 Refrigerated Microcentrifuge w/ COMBI-rotor	\$5,780