

Eppendorf Advantage

Advance your research through insights,
solutions, optimal outcomes, and
special offers.

JANUARY 1 – JUNE 30, 2026



Eppendorf Advantage Promotions

Advantage means putting yourself in a better position, and that's exactly what we want for you.

The work you do is more important now than ever before, and the samples you produce more precious. You have been able to rely on Eppendorf to help deliver accurate, reliable results, but we want to be more than your provider of consumables, instruments, and services.

We are also here to be your partner and resource center. Our bi-annual publication includes articles and promotional offers. We hope to help give you the Advantage to succeed in the lab!

Ordering Information

Product	Promo \$	Order no.
Centrifuges		
Save 38% on the Centrifuge 5427 R		
<i>Centrifuge 5427 R, Rotary knob with 30x1.5/2.0 mL rotor (FA-45-30-11)</i>	\$7,999	86-2965
Thermocyclers		
Save 38% on the Mastercycler® X50a		
<i>Mastercycler® X50a, aluminum block, 96-well plate or 0.1/0.2 mL tubes, 100-240 V/50-60 Hz</i>	\$7,995	86-2958
Save 38% on the Mastercycler® X40		
<i>Mastercycler® X40, aluminum block, 96-well plate or 0.1/0.2 mL tubes, 100-240 V/50-60 Hz</i>	\$4,995	86-2959
Incubator Shakers		
Save 35% on the New Brunswick S41i CO₂ Incubator Shaker with universal platform		
<i>New Brunswick S41i CO₂ incubator shaker, 1" orbit, 120 V/50-60 Hz, w/14"x 24" universal platform</i>	\$19,500	86-2960
Concentrators		
Save 39% on a Vacufuge® plus system with rotors for 2.0 mL, 15 mL, and 50 mL tubes		
<i>Vacufuge® plus system with integrated diaphragm vacuum pump, 120 V/50-60 Hz (US), with 48x1.5/2.0 mL rotor (F-45-48-11) & 6x15/50 mL rotor (F-35-6-30)</i>	\$9,999	86-2966
Mixers		
Get 33% off when you buy the ThermoMixer®C with ThermoTop®		
<i>ThermoMixer® C with ThermoTop®, 100-130 V/50-60 Hz</i>	\$4,500	86-2961
Pipettes		
Savings on the Eppendorf Research® 3 neo		
Get free color tag sample pack and ep Dualfilter T.I.P.S.® when you buy the new Eppendorf Research® 3 neo		
<i>Eppendorf Research® 3 neo 1-channel; 1 – 20 µL, yellow, variable volume; ACT-labeled. Includes 960 ep Dualfilter T.I.P.S.®, 2 – 20 µL, 53 mm, and ColorTag pipette marking rings sample bag.</i>	\$470	86-2962
<i>Eppendorf Research® 3 neo 1-channel; 10 – 200 µL, yellow, variable volume; ACT-labeled. Includes 960 ep Dualfilter T.I.P.S.®, 2 – 200 µL, 55 mm, and ColorTag pipette marking rings sample bag.</i>	\$470	86-2963
<i>Eppendorf Research® 3 neo 1-channel; 50 – 1000 µL, blue, variable volume; ACT-labeled. Includes 960 ep Dualfilter T.I.P.S.®, 50 – 1000 µL, 76 mm, and ColorTag pipette marking rings sample bag.</i>	\$470	86-2964

Eppendorf®, the Eppendorf Brand Design, ThermoMixer®, ThermoTop®, Mastercycler®, SmartBlock®, and Vacufuge® are registered trademarks of Eppendorf SE, Germany. Eppendorf New Lab Start Up Programs are subject to additional terms and conditions, view websites for full details. Offer valid in the U.S. only from January 1, 2026 through June 30, 2026 or while supplies last; void where prohibited by law or company policy. Product appearance, specifications, and/or prices are subject to change or cancellation at any time. All rights reserved including graphics and images. Copyright © 2026 by Eppendorf SE. ENA.C4.013.25