



Buy any two of the below products and
get a **FREE MiniMag Magnetic Stirrer**
with **promo code - BENCH MINIMAG**

Run Dates: 7/7/2025-9/30/2025

Buy two of these

Get 31-510 FREE



Cat #	Description of Product
31-190	hybex™ Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk.
31-191	hybex™ Media storage bottle, 100ml with standard (GL45) blue cap, 10/pk.
31-192	hybex™ Media storage bottle, 250ml with standard (GL45) blue cap, 10/pk.
31-193	hybex™ Media storage bottle, 500ml with standard (GL45) blue cap, 10/pk.
31-194	hybex™ Media storage bottle, 1000ml with standard (GL45) blue cap, 10/pk.
31-195	hybex™ Media storage bottle, 2000ml with standard (GL45) blue cap, 5/pk.
31-178	hybex™ Media storage bottle, 5L with standard (GL45) blue cap, 1/ea.
31-179	hybex™ Media storage bottle, 10L with standard (GL45) blue cap, 1/ea.
31-180	hybex™ Erlenmeyer Flask, 250ml with standard (GL45) blue cap, case of 12
31-181	hybex™ Erlenmeyer Flask, 500ml with standard (GL45) blue cap, case of 6
31-182	hybex™ Erlenmeyer Flask, 1000ml with standard (GL45) blue cap, case of 6
31-183	hybex™ Erlenmeyer Flask, 2000ml with standard (GL45) blue cap, 1 Each