



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## Liver Lysate

CATALOG NUMBER: 1404

## Specifications

|                     |        |
|---------------------|--------|
| Species             | Mouse  |
| Diagnosis           | Normal |
| Tested Applications | WB     |

## Properties

|                    |  |
|--------------------|--|
| Physical State     | Liquid   |
| Buffer             | Liver lysate is supplied in SDS sample buffer containing 50 mM DTT.  |
| Concentration      | 2 mg/mL  |
| Storage Conditions | Store at 2-8°C for continuous use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt. |

## Disclaimer

|            |   |
|------------|---|
| Disclaimer | This product is intended for laboratory research purposes only and should be used by qualified personnel only. ProSci lysates are not intended for use in humans. ProSci is not liable for damages or injuries resulting from receipt and/or use of ProSci lysates. |
|------------|---|

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).