

DIRAS3 Antibody


CATALOG NUMBER: 8139

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	DIRAS3 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human DIRAS3. The immunogen is located within the last 50 amino acids of DIRAS3.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	DIRAS3 antibody is human specific.
Predicted Molecular Weight	Predicted: 25 kDa Observed: 24 kDa

Properties

Purification	DIRAS3 antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	DIRAS3 antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	DIRAS3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



DIRAS3 Antibody 1

Western blot analysis of DIRAS3 in human testis tissue lysate with DIRAS3 antibody at 1 µg/ml.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).