

NALP1 Antibody


CATALOG NUMBER: 3037

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	NALP1 antibody was raised against a peptide corresponding to 13 amino acids near the carboxy terminus of human NALP1. The immunogen is located within the last 50 amino acids of NALP1.
Conjugate	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	NALP1 Antibody is immunoaffinity chromatography purified IgG.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	NALP1 Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	NALP1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.




NALP1 Antibody 1

Western blot analysis of NALP1 in U937 cell lysate with NALP1 antibody at (A) 1, (B) 2 and (C) 4 µg/mL.



NALP1 Antibody 2

Immunocytochemistry of NALP1 in K562 cells with NALP1 antibody at 10 µg/mL.



NALP1 Antibody 3

Immunofluorescence of NALP1 in K562 cells with NALP1 antibody at 20 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).