



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

XEDAR Antibody


CATALOG NUMBER: 3855

Specifications

Host Species	Rabbit
Species Reactivity	Human
Immunogen	XEDAR antibody was raised against recombinant human XEDAR.
Conjugate	Unconjugated
Tested Applications	ELISA, IHC-P, WB
User Note	Optimal dilutions for each application to be determined by the researcher.


Properties

Purification	XEDAR Antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	XEDAR Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	XEDAR antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



XEDAR Antibody 1

Western blot analysis of XEDAR in 293 cell lysate with XEDAR antibody at (A) 0.5, (B) 1 and (C) 2 µg/mL.



XEDAR Antibody 2

Immunohistochemistry of XEDAR in human skin tissue with XEDAR antibody at 10 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).