

Andes Virus Glycoprotein Antibody

CATALOG NUMBER: 7677

Specifications

Host Species	Rabbit
Species Reactivity	Virus
Immunogen	Andes virus glycoprotein antibody was raised against a 19 amino acid peptide near the amino terminus of the Andes virus glycoprotein. The immunogen is located within amino acids 70 - 120 of Andes Virus Glycoprotein.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.
Predicted Molecular Weight	N/A

Properties

Purification	Andes Virus Glycoprotein antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	Andes virus glycoprotein antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Andes virus glycoprotein antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).