

## SOX18 Antibody


CATALOG NUMBER: 7919

### Specifications

Host Species	Rabbit
Species Reactivity	Human, Mouse, Rat
Immunogen	SOX18 antibody was raised against a 17 amino acid peptide near the amino terminus of human SOX18. The immunogen is located within amino acids 40 - 90 of SOX18.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	SOX18 antibody is human, mouse and rat reactive. SOX18 is predicted to not cross-react with other SOX proteins.
Predicted Molecular Weight	Predicted: 42 kDa Observed: 44 kDa

### Properties

Purification	SOX18 antibody is affinity chromatography purified via peptide column.
Clonality	Polyclonal
Isotype	IgG
Physical State	Liquid
Buffer	SOX18 antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	SOX18 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



SOX18 Antibody 1

Western blot analysis of SOX18 in 3T3 cell lysate with SOX18 antibody at 1 µg/ml.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).