



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

ORAI3 Antibody [1B4F1]


CATALOG NUMBER: PM-4913

Specifications

Host Species	Mouse
Species Reactivity	Human, Rat
Immunogen	A 19 amino acid synthetic peptide from near the carboxy terminus of human ORAI3.
Conjugate	Unconjugated
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.

Properties

Purification	ORAI3 Monoclonal Antibody is immunoaffinity chromatography purified IgG.
Clonality	Monoclonal
Isotype	IgG2b
Physical State	Liquid
Buffer	ORAI3 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	ORAI3 monoclonal antibody can be stored at -20°C, stable for one year.



ORAI3 Antibody [1B4F1] 1

Western blot analysis of ORAI3 in rat lung tissue lysate with ORAI3 antibody at 2 µg/mL.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).