

DDDDK Antibody [1D1B12] (biotin)


CATALOG NUMBER: PM-7657-biotin

Specifications

Host Species	Mouse
Species Reactivity	Tag
Immunogen	Biotin-DDDDK antibody was raised against a synthetic peptide containing DYKDDDDK.
Conjugate	Biotin
Tested Applications	ELISA, WB
User Note	Optimal dilutions for each application to be determined by the researcher.
Specificity	Biotin-DDDDK-tag antibody will detect carboxy-terminal DDDDK-tagged recombinant proteins.
Predicted Molecular Weight	N/A

Properties

Purification	Biotin-DDDDK antibody is affinity chromatography purified via peptide column.
Clonality	Monoclonal
Isotype	IgG
Physical State	Liquid
Buffer	Biotin-DDDDK antibody is supplied in PBS containing 1% BSA and 0.02% sodium azide.
Concentration	1 mg/mL
Storage Conditions	Biotin-DDDDK antibody can be stored at 4°C for three months and -20°C, stable for up to one year.



DDDDK Antibody [1D1B12] (biotin) 1

Western blot analysis of (1) 1000ng, (2) 500ng, (3) 250ng, and (4) 125ng of DDDDK-tagged recombinant protein (GGP1) with Biotin-DDDDK-tag antibody at 0.5 µg/ml.

Disclaimer

Disclaimer	Optimal dilutions/concentrations should be determined by the end user. The information provided is a guideline for product use. This product is for research use only.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).