



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## TIM-3 ELISA Matched Antibody Pair (Risk Free)


CATALOG NUMBER: RFP-1001

### Specifications

Host Species	Mouse
Species Reactivity	Human
Immunogen	TIM-3 antibodies were raised against the extracellular domain of human TIM-3
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	Antibodies are supplied as protein A purified IgG1.
Clonality	Monoclonal
Isotype	IgG1
Physical State	Liquid
Buffer	PBS containing 0.02% sodium azide.
Concentration	Antibody 1 mg/mL
Storage Conditions	Antibodies can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.



**TIM-3 ELISA Matched Antibody Pair (Risk Free) 1**

A sandwich ELISA was performed using the anti-TIM-3 mAb RF16104 (5 µg/ml) as the capture antibody. Biotin-labeled anti-TIM-3 mAb RF16103-biotin (1 µg/ml) and streptavidin-HRP (0.1 µg/ml) were used...

**Disclaimer**

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).