



ProSci Incorporated  
12170 Flint Place  
Poway, CA 92064, USA  
prosci-inc.com

P: +1 (888) 513-9525  
Local: +1 (858) 513-2638  
Fax: +1 (858) 513-2692

customercare@prosci-inc.com  
techsupport@prosci-inc.com

## PD1 Detection Set

CATALOG NUMBER: SD8600

### Specifications

Host Species	Llama
Species Reactivity	Human
Immunogen	PD-1 antibodies were raised against a recombinant protein corresponding to amino acids 25 - 167 of human PD-1.
Conjugate	Unconjugated
Tested Applications	ELISA
User Note	Optimal dilutions for each application to be determined by the researcher.

### Properties

Purification	PD-1 Antibodies are affinity chromatography purified via Nickel column. Antibodies are supplied as a His-tagged purified protein. They also contains a myc-tag for detection.
Clonality	Monoclonal
Isotype	sdab
Physical State	Liquid
Buffer	PD-1 Antibodies are supplied in PBS.
Concentration	Antibody 1 mg/mL
Storage Conditions	PD-1 antibodies should be stored in working aliquots at -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Disclaimer

Disclaimer	This product is for research use only.
------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).